


Experiment 14

IR-NMR Exercise



Preparation

Before beginning this experiment, you should have

1. studied 'Structure Determination: Mass and Infrared Spectroscopy' and
2. studied 'Structure Determination: Nuclear Magnetic Resonance'.

You may also wish to read Chapter 29 in J.W. Zubrick's 'The Organic Chem Lab Survival Manual: A Students Guide to Techniques' pp.201-222, and Chapter 30 in J.W. Zubrick's 'The Organic Chem Lab Survival Manual: A Students Guide to Techniques' pp.223-233.

Objectives

Throughout the course, you have been exposed to the technique of 'Infrared Spectroscopy'. This has allowed you to correctly predict or confirm the presence or absence of certain functional groups in organic compounds. The purpose of this experiment is to test your ability to interpret and correlate IR and ^1H -NMR spectra data. From this information you will be able to identify functional groups and arrangements of hydrogen atoms within a molecule. Students are encouraged to discuss their approach to interpreting the spectra of the unknowns with their instructor and academic expert.


Introduction to Infrared Spectroscopy: Theory and Practice

Electromagnetic radiation

See 'Infrared Spectra Analysis Review' of the CHEM360 Lab Manual

Infrared radiation

Remember, infrared radiation carries relatively low levels of energy (e.g. ~ 1 -10 kcal/mol) which, when absorbed, result in only bond vibrations - stretching, rotating, bending and scissoring (i.e. deformation).



stretching vibration

bending vibration

Every molecule, depending on its makeup, is capable of absorbing infrared photons and increasing the intensity of its molecular motions. Different functional groups within the molecule will absorb photons at different infrared wavelengths. Thus, when a spectroscopic wavelength scan is performed on an organic molecule certain λ will be absorbed while other λ will pass through. Once we have the infrared spectrum of a compound, the spectrum can be analyzed and compared with known infrared absorptions for particular functional groups (see Table 14.1 and Table 14.2).

When analyzing a spectrum, it is important to look at four different regions of the spectrum for the presence or absence of specific absorption peaks. **Note:** you are not required to analyze the fingerprint region.

Table 14.1 Four Regions of the Infrared Spectrum

Wavenumber cm^{-1}				
4000	3000	2000	1400	600
N-H O-H	CH	C \equiv N C \equiv C	C=C C=O C=N	fingerprint region

Table 14.2 Correlation table of infrared absorption and functional group.

Type of Absorption	Wavenumber (cm ⁻¹)	Intensity of Absorption	Absorption of:
O-H stretch	3400-3640 2500-3300	strong, broad strong, very broad	alcohol carboxylic acid
N-H stretch	3310-3350	medium ('W' shape)	amine (1°)
C-H stretch	3300 3030 3020-3100 2850-2960 2750 & 2850	strong medium medium medium to strong weak-medium ('W' shape)	sp C-H of alkyne aromatic sp ² C-H of alkene sp ³ C-H of alkane O=C-H of aldehyde
C≡N stretch	2210-2260	medium, sharp	nitrile
C≡C stretch	2100-2260	medium, sharp	alkyne
C=O stretch	1670-1780 1730-1750 1720-1740 1705-1725 1700-1725 1640-1700 ca 1800 and 1760	strong, sharp	carbonyl ester aldehyde ketone carboxylic acid amide anhydride
C=C stretch	1650-1670 1600, 1500, 1450	weak-medium, sharp strong sharp	alkene aromatic
C=N stretch	1640-1670	medium, sharp	imine
N-H bend	1500-1650	medium to strong, sharp	amine and amide
N=O stretch	1500-1600 (1540) and 1320-1390	strong, sharp	nitro-compound
C-N stretch	1030, 1230	medium	amine
C-O stretch	1050-1150 1250-1310 1240 1175	strong strong broad strong, broad strong, broad	alcohol ester-conjugated ester-acetates ester-unconjugated
C-Cl stretch (terminal)	600-800	strong	alkyl halide
Ar-Cl stretch	1000-1175	medium-strong	aryl halide
C-Br stretch (terminal)	500-760	strong	alkyl halide
C-I (terminal)	500	strong	alkyl halide

Note: when a C=C bond is in conjugation with a carbonyl, the observed carbonyl absorption frequency will be $< \sim 30 \text{ cm}^{-1}$.

Calculation of the # degrees of unsaturation in a compound

(*See also 'Alkenes: Structure and Reactivity, Calculating a Molecules Degree of Unsaturation' McMurry 5th ed., p.190-192.

$$\text{Number of Degrees of Unsaturation} = nC + 1 + 1/2N - 1/2 nH - 1/2 nX$$

e.g., Therefore, for Compound A, $C_7H_{12} = (7) + 1 + 1/2(0) - 1/2 (12) - 1/2(0)$
 $= 7 + 1 - 6 = 2$ degrees of

unsaturation in Compound A. Note: an aromatic ring = 4 degrees of unsaturation, 1 for the ring + 3 for the 3 double bonds = 4

In this exercise, you will be provided with infrared spectra, and hints to the λ of interest.

NMR Spectroscopy—Theory and Practice

In physics courses you would learn about many properties of electrons including their mass, charge and spin. These properties may appear to be self-evident, but really they are only conceptual models to help describe their complex behaviour. Science trivia experts probably know the following properties of the electron:

diameter = $\sim 10^{-12}$ cm

rest mass = $m_e = 9.109534 \times 10^{-31}$ kg or 0.5110041 MeV or 1/1837 of H nucleus charge = 1.60219×10^{-19} C

specific charge = $e/m_e = 1.7588047 \times 10^{11}$ C kg⁻¹ magnetic moment = 9.284832×10^{-24} J T⁻¹

spin = $\frac{1}{2}$

Note: The property of 'spin' is not exactly like the earth rotating on its axis or a merry-go-round. Rather it tells us what a particle looks like from different directions. A particle of spin = $1/2$ must go through two rotations to look the same! Remember 'spin' is a quantum number which explains the splitting of spectral lines and that the 'spin' of a particle can line up parallel or anti-parallel with the magnetic field of an atom.

Like the electron, many nuclei also have spin properties (e.g., ^1H , ^{13}C , ^2H , ^{14}N , ^{19}F , ^{31}P). Those that do are all with odd numbered masses or even numbered masses with odd atomic numbers. Those that do have spin properties are known as magnetic nuclei. Thus, nuclei that have even masses and atomic numbers (e.g., ^{12}C , ^{16}O , ^{32}S) are non-magnetic.

Magnetic nuclei will have nuclear spins without a specific orientation, if not in a strong external magnetic field. However, when placed into a magnetic field the nuclei will orient themselves parallel (favoured lower energy state) or anti-parallel (less favoured higher energy state) with the field. Then if the nuclei are exposed to radio waves of the right frequency, energy absorption can occur, and the parallel spin will convert to anti-parallel spin (spin-flip). This spin-

flipping between energy states of the two spin orientations of the nucleus is what is meant by 'nuclear magnetic resonance'.

Nuclear magnetic resonance (NMR) spectroscopy is a very useful tool for organic chemists. Used in conjunction with mass and IR spectroscopy, it allows you to form a framework of the carbon-hydrogen and carbon-carbon bonds in a molecule from which you can infer the structure and identity of a compound.

When you look at a ^1H -NMR spectrum you are trying to decipher and correlate four different kinds of information. Because the information in NMR spectra sometimes is very complex, you must develop a systematic approach to analyzing NMR spectra.

1. First the spectrum contains **chemical shift (δ) information** (i.e. the position of the peak) which tells you about the structural grouping to which the H is bound (analogous to relationship between infrared λ and functional groups). Chemical shifts (δ) are measured in **hertz (Hz) or ppm** and are the distance from the center of the signal to a reference signal, usually tetramethylsilane (TMS). Refer to the chemical shifts shown in Table 14.3.
2. Second, the size of or area under a peak (determined by **integration**) tells you about number of identical H in a particular electronic environment.

Table 14.3 ^1H NMR Chemical Shifts (σ , ppm) for Various Functional Groups

σ , ppm			σ , ppm	
TMS $(\text{CH}_3)_4\text{Si}$	0	TMS is used as reference for both ^1H and ^{13}C -NMR. It gives rise to a single peak that occurs upfield (farther right) of other absorptions normally found in organic compds.	Alcohols, ethers	
Cyclopropane	0.0-0.4		$\begin{array}{c} \\ \text{HO}-\text{C}-\text{H} \\ \end{array}$	3.3-4.0
Alkanes RCH_3 R_2CH_2 R_3CH	0.7-1.3 1.2-1.4 1.4-1.7		$\begin{array}{c} \\ \text{RO}-\text{C}-\text{H} \\ \end{array}$	3.3-4.0
Alkenes $\begin{array}{c} \quad \\ -\text{C}=\text{C}-\text{H} \text{ (vinyl)} \\ \quad \end{array}$ $\begin{array}{c} \quad \\ -\text{C}=\text{C}-\text{CH}_3 \text{ (allyl)} \\ \quad \end{array}$	4.6-5.9 1.6-1.9	Vinylic protons are strongly deshielded by the neighbouring pi bond and therefore absorb in this characteristic downfield position. protons on C next to unsaturated (allylic, benzylic, next to carbonyl) show charact. absorptions in this region, just downfield from other alkane resonance.	Esters $\begin{array}{c} \text{O} \\ \\ \text{R}-\text{C}-\text{O}-\text{C}-\text{H} \\ \end{array}$ $\begin{array}{c} \text{O} \\ \\ \text{RO}-\text{C}-\text{C}-\text{H} \\ \end{array}$	3.7-4.1 2.0-2.6
Alkynes $-\text{C}\equiv\text{C}-\text{H}$ (alkynyl) $-\text{C}\equiv\text{C}-\text{CH}_3$	2.5-2.7 1.8		Carboxylic acids $\begin{array}{c} \text{O} \\ \\ \text{HO}-\text{C}-\text{C}-\text{H} \\ \end{array}$ $\begin{array}{c} \text{O} \\ \\ \text{R}-\text{C}-\text{O}-\text{H} \end{array}$	2.0-2.6 10.5-12
Aromatic $\text{Ar}-\text{H}$ (aryl) $\begin{array}{c} \\ \text{Ar}-\text{C}-\text{H} \text{ (benzyl)} \\ \end{array}$	6.5-8.0 2.5-2.7	protons on aromatic rings (aryl protons) are strongly deshielded by the pi orbitals of the ring and show charact. absorptions in this lower-field range.	Ketones $\begin{array}{c} \text{O} \\ \\ \text{R}-\text{C}-\text{C}-\text{H} \\ \end{array}$	2.1-2.4
Fluorides, $\begin{array}{c} \\ \text{F}-\text{C}-\text{H} \\ \end{array}$	4.0-4.5	protons on C next to electronegative atoms (X, O, N) are deshielded because of the electron-withdrawing ability of these atoms. Thus the protons absorb in this midfield range.	Aldehydes $\begin{array}{c} \text{O} \\ \\ \text{R}-\text{C}-\text{H} \end{array}$	9.7-10.0
Chlorides, $\begin{array}{c} \\ \text{Cl}-\text{C}-\text{H} \\ \end{array}$	3.0-4.0		Amides $\begin{array}{c} \text{O} \\ \\ \text{R}-\text{C}-\text{N}-\text{H} \\ \end{array}$	5.0-8.0
$\begin{array}{c} \text{Cl} \\ \\ \text{Cl}-\text{C}-\text{H} \\ \end{array}$	5.8		Alcohols $\text{R}-\text{O}-\text{H}$ 2.5-5.0 Phenols $\text{Ar}-\text{O}-\text{H}$ 4.0-12.0 Amines $\text{R}-\text{NH}_2$ 1.0-5.0	
Bromides, $\begin{array}{c} \\ \text{Br}-\text{C}-\text{H} \\ \end{array}$	2.5-4.0	protons on oxygen bearing C atoms are deshielded by the electron-withdrawing effect of the nearby O. Splitting of OH proton not usually observed, therefore usually seen as a broad singlet.		
Iodides, $\begin{array}{c} \\ \text{I}-\text{C}-\text{H} \\ \end{array}$	2.0-4.0			
Nitroalkanes, $\begin{array}{c} \\ \text{O}_2\text{N}-\text{C}-\text{H} \\ \end{array}$	4.2-5.6			


The integration data may be in the form of 'peak areas' or an 'integration line' on the NMR spectrum. If the integration data is in the form of peak areas, the number of protons responsible for signal is easily obtained from the integration by looking at the molecular formula for the unknown and using the following equation:

$$\# \text{ of Protons Responsible for Signal} = \text{total \# of protons in molecule} \times \frac{\text{area under signal}}{\sum \text{ area under all}}$$

signals

e.g., If the molecular formula of the unknown is C_6H_{11} , and has 3 signals with areas of 31, 6, and 20 (sum = 57), then the number of protons responsible for the first signal is $11 \times 31/57 = 5.98 \approx 6$; and $11 \times 6/57 = 1.15 \approx 1$ for the second signal; and $11 \times 20/57 = 3.86 \approx 4$ for the third signal ($6 + 1 + 4 = 11$).

If the integration data is in the form of an 'integration line' on the NMR spectrum, the number of protons responsible for signal is obtained by measuring the height (h) of each of the integration lines for each of the signals. The **ratio of the signal's heights** is then compared with the molecular formula to determine the number of protons responsible for each signal.



Note: in this experiment, the integration data has been analyzed and given to you.

- The third type of information you are looking for is called **signal splitting or multiplicity**. This information tells you about the number of neighbouring H and is calculated using the 'n+1 Rule'.

$$N = n + 1$$

where N = number of peaks observed, and n = number of equivalent adjacent H (see Table 14.4 for calculating signal splitting).

Table 14.4 Signal Splitting Calculation Using Pascal's Triangle ($N=n+1$)

N = number of peaks observed for absorbing protons.

n = number of equivalent adjacent hydrogens.

	<u>n+1</u>	<u>n</u>	splitting
1	1	0	singlet
1 1	2	1	doublet
1 2 1	3	2	triplet
1 3 3 1	4	3	quartet
1 4 6 4 1	5	4	quintet
1 5 10 10 5 1	6	5	sextet
1 6 15 20 15 6 1	7	6	septet

Pascal's Triangle

- The last type of information to seek has to do with **coupling constants (J)**.

This is the distance between two adjacent peaks in the signal and is measured in hertz. For ideal triplets and quartets, their signal peaks should be symmetrical with relative peak areas of 1:2:1 and 1:3:3:1 respectively. In actual spectra this is rarely the case, and the peak ratios are distorted.

In Figure 14.5, the peaks of two sets of protons are interacting with each other, i.e., they share the **same** coupling constant (J_{ab}), they are leaning towards each other, and the peaks on the side nearest the other set of signals are higher than predicted.




Figure 14.5 Coupling constant (J) for nearest neighbouring H's.

The protons of 'set H_a ' have split the protons of 'set H_b ' into a quartet and the protons of 'set H_b ' have split the protons of 'set H_a ' into a triplet. The interacting triplet-quartet peaks of this kind are evidence of the presence of an ethyl group (CH_3CH_2-) in the molecule.

In summary, you must systematically analyze a ^1H -NMR and correlate four different kinds of information: the chemical shift, the integration, the multiplicity or signal splitting pattern, and the coupling constants.

Finally, as a chemistry student, it is essential that you have a good understanding of the concepts and practice of ^1H -NMR. Unfortunately, modern NMR spectrometers are very expensive to own and require a great deal of technical experience to operate. Therefore, in this course you will only be able to gain experience in the theory, interpretation and analyses of ^1H -NMR spectra. This exercise, where you are provided with unknown ^1H -NMR spectra plus integration data, will allow you learn and practice this important area of chemistry.

Procedure

1. Do the instructor led and practice problems.
2. Obtain 4 unknown samples of organic compounds, and perform melting points or boiling points, and mass, infrared, and $^1\text{H-NMR}$ spectroscopy on each of them (This may have been already done for you!).
3. Analyze the provided spectral data and determine the identity and structure of each unknown. (Your instructor will provide you with a handout for each unknown).
4. Calculate the 'Degrees of Unsaturation' present in the molecule.
Number of Degrees of Unsaturation = $n\text{C} + 1 + 1/2n\text{N} - 1/2 n\text{H} - 1/2 n\text{X}$
5. Then analyze the infrared spectrum.
6. It is suggested you analyze your $^1\text{H-NMR}$ spectra as follows:
 - a) Always survey the spectrum from left (downfield) to right (upfield). As a rule, the further the H signal is shifted downfield, the more electronegative the group to which the H is bound.
 - b) Number/code the signals, furthest downfield being #1 or A.
 - c) Ask yourself how many hydrogens are giving rise to each signal?
 - d) Determine the chemical shift of the signal (δ), in PPM.
 - e) Next determine the multiplicity of the signal. Is it a singlet, doublet, triplet, quartet, sextet, septet, multiplet? From this you will be able to say how many H are on the neighbouring carbon.
 - f) Now begin to build the framework of C, H, O, and N, which is also supported by the information obtained from the mass and infrared spectra. From this you will eventually obtain the complete structure of the molecule.
 - g) Once you have arrived at a structure, confirm the presence (or absence) or the expected peaks such a molecule would give in the Infrared and $^1\text{H-NMR}$ spectra. Check to see if the other physical data matches your conclusion (mp or bp, and chemical formula).
7. Submit your answers, along with your reasoning and spectra, in a brief lab report (see 'Write-up Instructions'). If you have any questions about the data, please call your instructor or AE.

How to Interpret an Infrared Spectrum

Step 1 Divide the infrared spectrum into four main areas (use pencil and ruler and take into account any off shift in the spectrum's wavenumbers).

- i) Above 3000 cm^{-1}
- ii) Between 3000 and 2000 cm^{-1}
- iii) Between 2000 and 1400 cm^{-1}
- iv) Below 1400 cm^{-1} (fingerprint region)

Step 2 Starting at the left of the spectrum, examine the area **above 3000 cm^{-1}** , first looking in the region near 3300 cm^{-1} and record in tabular format the presence/absence of:

- i) A broad, very strong absorption band of an '**O-H**'. If present, it means you know that your molecule is at least an **alcohol**.
- ii) A broad, weak to medium strength, double or single absorption band of '**N-H**'. If present it means you have an **amine** (1° or 2°) or possibly an **amide**.
- iii) A sharp, medium to strong, single absorption band of ' **$\equiv\text{C-H}$** ' of a **terminal alkyne**. Note: If present, it means you should also see a ' **$\text{C}\equiv\text{C}$** ' absorption near 2250 cm^{-1} .

After examining the region around 3300 cm^{-1} , look for any sharp, weak to medium absorption just above 3000 cm^{-1} (e.g. 3050 cm^{-1}) resulting from the '**C-H**' stretch of a sp^2 hybridized carbon. If present, it means you have a ' **$\text{C}=\text{C-H}$** ' of an alkene or aromatic compound.

Step 3 Next examine the area between **3000 and 2000 cm^{-1}** and record the presence/absence of absorption bands or peaks.

- i) First look just below 3000 cm^{-1} (e.g. $2850\text{-}2950\text{ cm}^{-1}$) resulting from the '**C-H**' stretch of a sp^3 hybridized carbon. If present, it means you are seeing the '**C-H**' stretch of an **$-\text{CH}_2$ or $-\text{CH}_3$** group. Note: This absorption is not very informative as most organic compounds have $-\text{CH}_2$ or $-\text{CH}_3$ groups.
- ii) Then look for the extremely broad peak, starting at 3300 cm^{-1} and extending all the way to $\sim 2500\text{ cm}^{-1}$, caused by the **O-H dimer** between two **carboxylic acid** molecules (COOH). This absorption is probably the most difficult to see as other absorption peaks may be overlapping the broad peak.
- iii) Finally look for a sharp, weak to medium peak caused by either ' **$\text{C}\equiv\text{C}$** ' or ' **$\text{C}\equiv\text{N}$** '.
- iv) If present, then the compound is an alkyne (might also have the '**C-H**' of a terminal alkyne, see step 2 above) or a nitrile.

Step 4 Next examine the area between **2000 and 1400 cm^{-1}** and record the presence/absence of absorption bands or peaks.

- i) First look near 1700 cm^{-1} (e.g. $1680\text{-}1750\text{ cm}^{-1}$) for a sharp, strong peak resulting from the ' **$\text{C}=\text{O}$** ' stretch of a **carbonyl**. Note: This absorption is very informative and will be present if your compound is an aldehyde, ketone, ester, amide, or carboxylic acid.
- ii) Next look near 1650 cm^{-1} (e.g. $1600\text{-}1670\text{ cm}^{-1}$) for a sharp, weak peak resulting from the ' **$\text{C}=\text{C}$** ' stretch of an **alkene**.

- iii) Finally look near 1600 cm^{-1} and 1500 cm^{-1} for a sharp, double peak resulting from the '**C=C**' stretch of an **aromatic ring**.

Step 5 If you dare, you may look in the **fingerprint region (area below 1400 cm^{-1})** and record the presence of absorption bands or peaks.

- i) First look near 1200 cm^{-1} for a sharp, strong peak resulting from the '**C-O**' stretch of an **ester**. Note: This absorption is very difficult to see and may or may not be present, i.e. conclusive if present, inconclusive if not present.
- ii) If you suspect you have an aromatic ring (absorption bands at ~ 3030 and 1600 and 1500 cm^{-1} present), you may try to discern the substitution pattern of the benzene ring by looking at the strong absorption bands of the **ring 'C-H'** out-of-plane bending vibrations in the region $680-900\text{ cm}^{-1}$.

Benzene Substitution Pattern	Ring 'C-H' Absorption Bands Present (cm^{-1})
monosubstituted	2 sharp peaks, 730-770, 690-710
<i>ortho</i> disubstituted	1 sharp peak, 735-770
<i>meta</i> disubstituted	3 sharp peaks, 860-900, 750-810, 680-725
<i>para</i> disubstituted	1 sharp peak, 800-860
1,2,3 trisubstituted	2 sharp peaks, 760-780, 705-745
1,3,5 trisubstituted	2 sharp peaks, 810-865, 675-730
1,2,4 trisubstituted	2 sharp peaks, 870-885, 805-825

Ref: McMurry, J., 2000. Organic Chemistry, 5th ed, Brooks/Cole, p.578-579, (4th ed, p.559) Nakanishi, K., 1964. Infrared Absorption Spectroscopy, Holden Day p.27.

- iii) If you suspect you have a long straight chain ($>4\text{ C}$) alkane, (absorption bands at $28050-2950\text{ cm}^{-1}$ present but not much else), you may try to see the sharp, weak absorption due to the concerted rocking or $>4\text{-CH}_2$ in a chain. It lies in the region $720 \pm 10\text{ cm}^{-1}$.

Step 6 Finally, you will summarize your results by making a statement as to what functional groups you suspect to be present in the molecule or perhaps you will be asked to select from a list of suggested structures, which molecule most likely would generate the spectrum just analyzed.

As tabular format that you might find useful for recording your findings is shown below:

Infrared Data:

	Absorption Band#	Frequency (cm^{-1})	Peak Shape (sharp, broad)	Peak Intensity (strong, med. or weak)	Functional Group Indicated
$> 3000\text{ cm}^{-1}$					
Between 3000 and 2000 cm^{-1}					
Between 2000 and 1400 cm^{-1}					
$< 1400\text{ cm}^{-1}$					

Functional Group(s) absent:

How to Interpret a ^1H -NMR Spectrum

- Step 1** Survey the ^1H -NMR spectrum from left (downfield) to right (upfield). As a rule, the further the H signal is shifted downfield, the more electronegative the group to which the H is bound. Assign a code # for each signal or group of signals, furthest downfield being #1. Remember, each signal or group of signals represents a proton(s) in a different chemical environment.
- Step 2** Determine the **chemical shift** (δ or delta scale, PPM) for each signal, or group of signals, by measuring the distance from the reference TMS peak to the center of the signal or group of signals.
- Step 3** Determine the **number of hydrogens** giving rise to each signal or group of signals. This information may be provided (integration data), or you may have to measure the area under the peak(s) of each signal. It is important to check that the sum of the # of protons from each signal adds up to the total number of H indicated in the molecular formula.
- Step 4** Determine the **multiplicity or splitting pattern** of the signal or group of signals by counting the number of distinguishable peaks. This may be difficult at times as very weak peaks may be overshadowed by background noise. Also note that you can have irregular signals, and signals that have been split by several other nonequivalent protons. In this case, more advanced methods of analysis than shown here are then required to determine the multiplicity.
- Step 5** Use the **$N = n + 1$** rule to determine the number of **neighbouring hydrogens** for each signal. (N = number of peaks observed, and n = number of equivalent adjacent H.)

Splitting Pattern	# of Neighbouring Hydrogens
singlet	0
doublet	1
triplet	2
quartet	3
quintet	4
sextet	5
septet	6

see Table 14.4 for calculating signal splitting

- Step 6** Make your **signal assignments** after consulting the Table 14.3 ^1H -NMR Shifts for Various Functional Groups.
- Step 7** This is the most crucial step in the analysis, and the most difficult. Now you must begin to put together all your results, using the signal assignments, the infrared spectra data, the degrees of unsaturation, and molecular formula to piece together the structure of the unknown organic compound. The tabular format that you will find useful for recording your results is shown below:

Table #X. ^1H -NMR Spectral Data:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment

In the following pages, we present sample spectrum interpretations for 2-propanol and for 1-propanol.

Sample Interpretation of a ¹H-NMR Spectrum


Use the following examples to understand the step-by-step approach of how to use a ¹H-NMR spectrum to identify an unknown. Please note that the 'degrees of unsaturation' has been calculated, and the infrared spectrum has already been analyzed.

Experiment 14 Student Name: **Earle L. Mier** AU ID#: **2550100**

Unknown #: 360-14-A Mol. Wt.: **60.1 g/mol** Formula: **C₃H₈O** BP: **82.4° C**

Infrared Spectrum of 360-14-B (thin film)
CDCl₃)

¹H-NMR Spectrum of 360-14-B (90 MHz in



$$\begin{aligned} \text{Number of Degrees of Unsaturation} &= nC + 1 + 1/2nN \\ - 1/2 nH - 1/2 nX & \\ \text{for Compound 360-14-A, C}_3\text{H}_8\text{O} &= (3) + 1 + 1/2(0) \\ - 1/2 (10) - 1/2(0) & \\ &= (3 + 1 + 0 - 4 - 0) = \mathbf{0} \text{ degrees of unsaturation} \end{aligned}$$


Infrared Spectrum Data for Unknown 360-14-A:

	Absorption Band	Frequency (cm ⁻¹)	Peak Shape (sharp, broad)	Peak Intensity (strong, med. or weak)	Functional Group Indicated
> 3000 cm ⁻¹	A	3346	Broad	Strong	O-H stretch of alcohol
Between 3000 and 2000 cm ⁻¹	B	2972, 2933, 2884	Multiple, Sharp	Strong	sp ³ C-H stretches of alkane
Between 2000 and 1400 cm ⁻¹					
< 1400 cm ⁻¹					

Functional Group(s) absent: no N-H, C≡C-H, sp² C-H, C≡C, C≡N, C=O, C=C alkene or aromatic

Interpretation of the ¹H-NMR Spectrum for Unknown 360-14-A


Step 1 Survey the ¹H-NMR spectrum from left (downfield) to right (upfield), assigning a code # for each signal or group of signals, furthest downfield being #1. Two proton signals are shifted downfield, indicating they are 'close' to an electronegative group (deshielded). The spectrum also indicates that there are proton(s) in a total of 3 different chemical environment. Note that Signal 1 has been expanded for easier viewing.



¹H-NMR Spectral Data for Unknown 360-14-A:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1				Shifted downfield		
2				Shifted downfield		
3						


Step 2 Determine the **chemical shift** (δ or delta scale, PPM) for each signal, or group of signals, by measuring the distance from the reference TMS peak (0.0) to the center of the signal or group of signals.



¹H-NMR Spectral Data for Unknown 360-14-A:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	4.0			Shifted downfield		
2	2.15			Shifted downfield		
3	1.2					


Step 3 The **number of hydrogens** giving rise to each signal or group of signals may be provided as integration data, or you may have to determine the ratio of the area under the peak(s) of each signal. Use a ruler to measure the height of each integration line (from left to right = 0.5 cm : 0.6 cm : 2.8 cm). The sum of the # of protons from all signals must add up 8 = the total number of H indicated in the molecular formula. In this case, this is accomplished by multiplying by 2 and rounding off the ratio of heights (from left to right = 1 : 1 : 6) = 8 H total.



¹H-NMR Spectral Data for Unknown 360-14-A:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	4.0	1		Shifted downfield		
2	2.15	1		Shifted downfield		
3	1.2	6				

Step 4 Determine the **multiplicity or splitting pattern** of the signal or group of signals by counting the number of distinguishable peaks in each signal. Signal 1 is a septet (as shown on the expanded area below). Signal 2 is a singlet. Signal 3 is a doublet.



¹H-NMR Spectral Data for Unknown 360-14-A:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	4.0	1	septet	Shifted downfield		
2	2.15	1	singlet	Shifted downfield		
3	1.2	6	doublet			

Step 5 Use the $N = n + 1$ rule to determine the number of **neighbouring hydrogens** for each signal. (N = number of peaks observed, and n = number of equivalent adjacent H.)
 For Signal 1, $N = n + 1 = (7 = n + 1)$, so $n = 6$. Signals 2-4 # Neighbouring H are calculated the same way.

¹H-NMR Spectral Data for Unknown 360-14-A:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	4.0	1	septet	Shifted downfield	6	
2	2.15	1	singlet	Shifted downfield	0	
3	1.2	6	doublet		1	

Step 6 Make your **signal assignments**, consulting Table 14.3 ¹H-NMR Shifts for Various Functional Groups.

-According to Table 14.3, **Signal 1** (downfield 1H septet at δ 4.0) must be an alkyl halide, an ether, or an alcohol. It can't be an alkyl halide, because there is no halogen in the molecular formula, so it must be an ether, or an alcohol. It is therefore given the signal assignment of R-O-CH-R₂.

-We know that there is an alcohol in this unknown compound from the infrared analysis (absorption band A). Therefore, we must look to find a decoupled (0 neighbours), deshielded (shifted downfield) signal due to a single proton bonded directly to an oxygen. Only **Signal 2** (singlet, 1H at δ 2.15) could be the H of the hydroxyl group of an alcohol. It is therefore given the signal assignment of R-O-H.

-**Signal 3** (6H at δ 1.2) is due to an alkane R-CH₃ (δ 0.7-1.3), or R₂-CH₂ (δ 1.2-1.4), but not R₃-CH (δ 1.4-1.7). It is likely that Signal 3 is slightly shifted down field, and so it is more likely the signal assignment should be R-CH₃.

¹H-NMR Spectral Data for Unknown 360-14-A:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	4.0	1	septet	Shifted downfield	6	R-O-CH-R ₂
2	2.15	1	singlet	Shifted downfield	0	R-O-H of an alcohol
3	1.2	6	doublet		1	R-(CH ₃) ₂


Your ¹H-NMR Spectral Data Table is now complete. Proceed to Step 7 on the next page. Please note that you have not been asked to extract information from the spectrum regarding coupling constants (J).

Step 7 Again, this is the crucial step of the analysis. Let's begin to put together all your results by using the signal assignments, the infrared spectra data (O-H in molecule), the degrees of unsaturation (0), and molecular formula (C₃H₈O).

Rationale:


It was immediately apparent that unknown 360-14-A is an **alcohol**, -OH, due to the broad infrared absorption band O-H at 3346 cm⁻¹, and an downfield lone proton (Signal 2 singlet at δ 2.15).

(Note: 'Atomic Accounting' is now: C₃H₈O - OH = C₃H₇).




Signal 1 indicates there is a R-O-CH- (δ = 4.0, Integration = 1H, splitting pattern = septet) in the molecule. Also this O-CH- must have a neighbouring carbons with 6 H on it, so it must be O-CH- (CH₃)₂. This is confirmed by the presence of a R-CH₃ group (δ =1.2, Integration = 6H, splitting pattern = doublet). The last piece of the puzzle is that the CH₂ group is shifted downfield (δ =3.67) and so must be bonded to the only oxygen in the molecule, i.e. an O-CH-(CH₃)₂ group.

(Note: 'Atomic Accounting' is now: C₃H₇ - CH = C₂H₆).




Signal 3 confirms that there are 2 **methyl groups** R-CH₃ ($\delta=1.2$, Integration = 6H, splitting pattern = doublet). Also, these 2 -CH₃- both have the same neighbouring carbon with a single H, so Signal 3 must be due to (CH₃)₂-CH-.
 (Note: 'Atomic Accounting' is now: C₂H₆ - 2(CH₃) = 0, all atoms accounted for!).



In summary unknown #360-14-A contains an alcohol, and 2 methyl groups. Therefore, the **Unknown #360-14-A** must be **2-propanol**, which also has the correct/matching molecular weight/formula, and boiling point. Also, the ¹H-NMR Spectrum shows there are 3 different types of H environments present in the molecule (and all 3 have been accounted for).

Note: The calculation for the 'degrees of unsaturation' in this problem also makes sense. There are 0 degrees of unsaturation present in the unknown.



Another Sample Interpretation of a ¹H-NMR Spectrum


Use the following examples to understand the step-by-step approach of how to use a ¹H-NMR spectrum to identify an unknown. Please note that the 'degrees of unsaturation' has been calculated, and the infrared spectrum has already been analyzed.

Experiment 14 Student Name: R.B. Flask AU ID#: 2550100

Unknown #: 360-14-B Mol. Wt.: 60.1 g/mol Formula: C₃H₈O BP: 97° C

Infrared Spectrum of 360-14-B (thin film)
CDCl₃)

¹H-NMR Spectrum of 360-14-B (90 MHz in



$$\begin{aligned}
 \text{Number of Degrees of Unsaturation} &= nC + 1 + 1/2nN - 1/2 nH - 1/2 nX \\
 \text{for Compound 360-14-B, C}_3\text{H}_8\text{O} &= (3) + 1 + 1/2(0) - 1/2 (10) - 1/2(0) \\
 &= (3 + 1 + 0 - 4 - 0) = 0 \text{ degrees of unsaturation}
 \end{aligned}$$


Infrared Spectrum Data for Unknown 360-14-B:

	Absorption Band	Frequency (cm ⁻¹)	Peak Shape (sharp, broad)	Peak Intensity (strong, med. or weak)	Functional Group Indicated
> 3000 cm ⁻¹	A	3333	Broad	Strong	O-H stretch of alcohol
Between 3000 and 2000 cm ⁻¹	B	2963, 2938, 2878	Multiple, Sharp	Strong	sp ³ C-H stretches of alkane
Between 2000 and 1400 cm ⁻¹					
< 1400 cm ⁻¹					

Functional Group(s) absent: no N-H, C≡C-H, sp² C-H, C≡C, C≡N, C=O, C=C alkene or aromatic

Interpretation of the $^1\text{H-NMR}$ Spectrum for Unknown 360-14-B


Step 1 Survey the $^1\text{H-NMR}$ spectrum from left (downfield) to right (upfield), assigning a code # for each signal or group of signals, furthest downfield being #1. Two proton signals are shifted downfield, indicating they are 'close' to an electronegative group (deshielded). The spectrum also indicates that there are proton(s) in a total of 4 different chemical environment.



$^1\text{H-NMR}$ Spectral Data for Unknown 360-14-B:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1				Shifted downfield		
2				Shifted downfield		
3						
4						


Step 2 Determine the **chemical shift** (δ or delta scale, PPM) for each signal, or group of signals, by measuring the distance from the reference TMS peak (0.0) to the center of the signal or group of signals.



¹H-NMR Spectral Data for Unknown 360-14-B:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	3.6			Shifted downfield		
2	2.25			Shifted downfield		
3	1.6					
4	0.95					

Step 3 The number of hydrogens giving rise to each signal or group of signals may be provided as integration data, or you may have to determine the ratio of the area under the peak(s) of each signal. Use a ruler to measure the height of each integration line (from left to right = 1cm: 0.5cm: 1cm: 1.5cm). The sum of the # of protons from all signals must add up 8 = the total number of H indicated in the molecular formula. This is accomplished by multiplying the ratio of heights by 2 (from left to right = 2 : 1 : 2 : 3) = 8 H total.




1 2 3

¹H-NMR Spectral Data for Unknown 360-14-B:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	3.6	2		Shifted downfield		
2	2.25	1		Shifted downfield		
3	1.6	2				
4	0.95	3				

Step 4 Determine the **multiplicity or splitting pattern** of the signal or group of signals by counting the number of distinguishable peaks in each signal. Signal 1 is a triplet (as shown on the spectrum below). Signal 2 is a singlet. Signal 3 is difficult to interpret, but it is a quintet. Signal 4 is clearly a triplet.



¹H-NMR Spectral Data for Unknown 360-14-B:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	3.6	2	triplet	Shifted downfield		
2	2.25	1	singlet	Shifted downfield		
3	1.6	2	quintet			
4	0.95	3	triplet			

Step 5 Use the $N = n + 1$ rule to determine the number of **neighbouring hydrogens** for each signal. (N = number of peaks observed, and n = number of equivalent adjacent H.) For Signal 1, $N = n + 1 = (3 = n + 1)$, so $n = 2$. Signals 2-4 # Neighbouring H are calculated the same way.

¹H-NMR Spectral Data for Unknown 360-14-B:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	3.6	2	triplet	Shifted downfield	2	
2	2.25	1	singlet	Shifted downfield	0	
3	1.6	2	quintet		4	
4	0.95	3	triplet		2	

Step 6 Make your signal assignments, consulting Table 14.3 ¹H-NMR Shifts for Various Functional Groups.

-According to Table 14.3, **Signal 1** (2H at δ 3.6) must be an alkyl halide, an ether, or an alcohol. It can't be an alkyl halide, because there is no halogen in the molecular formula, so it must be an ether, or an alcohol. It is therefore given the signal assignment of R-O-CH₂-.

-We know that there is an alcohol in this unknown compound from the infrared analysis (absorption band A). Therefore we must look to find a decoupled (0 neighbours), deshielded (shifted downfield) signal due to a single proton bonded directly to an oxygen. Only **Signal 2** (singlet, 1H at δ 2.25) could be the H of the hydroxyl group of an alcohol. It is therefore given the signal assignment of R-O-H.

-According to Table 14.3, **Signal 3** (2H at δ 1.6) must be due to an alkane or alkene. Since there is no C=C in the molecule (see infrared data, and degrees of unsaturation = 0), it must be due to an alkane R-CH₃, R₂-CH₂ or R₃-CH. Of the three possibilities, the table suggests it could be R₃-CH (δ 1.4-1.7). But what if Signal 3 is near enough to be coupled with protons bonded to a carbon bonded to an oxygen. They might be slightly deshielded as well, due their 'close' proximity to the electron- withdrawing oxygen. So this signal could also be due to R₂-CH₂ (δ 1.2-1.4). This makes more sense, since we are

dealing with a 2H signal (integration = 2). Thus Signal 3 is given the signal assignment R₂-CH₂.

-Signal 4 (3H at δ 0.95) must be due to an alkane R-CH₃ (δ 0.7-1.3).

¹H-NMR Spectral Data for Unknown 360-14-B:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	3.6	2	triplet	Shifted downfield	2	R-O-CH ₂ -
2	2.25	1	singlet	Shifted downfield	0	R-O-H of an alcohol
3	1.6	2	quintet		4	R ₂ -CH ₂
4	0.95	3	triplet		2	R-CH ₃

Your ¹H-NMR Spectral Data Table is now complete. Proceed to Step 7.


Please note that you have not been asked to extract information from the spectrum regarding coupling constants (J).

Step 7 This is the most crucial step in the analysis, and the most difficult. Now you must begin to put together all your results by using the signal assignments, the infrared spectra data (O-H in molecule), the degrees of unsaturation (0), and molecular formula (C₃H₈O).

Rationale:


It was immediately apparent that unknown 360-14-B is an **alcohol**, -OH, due to the broad infrared absorption band O-H at 3333 cm⁻¹, and a downfield lone proton (Signal 2 singlet at δ 2.25).

(Note: 'Atomic Accounting' is now: C₃H₈O - OH = C₃H₇).




Signal 1 indicates there is a R-O-CH₂- (δ = 3.6, Integration = 2H, splitting pattern = triplet) in the molecule. Also, this O-CH₂- must have a neighbouring carbon with 2 H on it, so it must be O-CH₂-CH₂. This is confirmed by the presence of another -CH₂ group (δ=1.6, Integration = 2H, splitting pattern = quintet). The last piece of the puzzle is that the CH₂ group is shifted downfield (δ=3.67) and so must be bonded to the only oxygen in the molecule, i.e. an **O-CH₂-CH₂- group**.

(Note: 'Atomic Accounting' is now: C₃H₇ - C₂H₄ = CH₃).




Signal 3 confirms that there is a methylene group R₂CH₂ (δ=1.6, Integration = 2H, splitting pattern = quintet). Also, this -CH₂- has neighbouring carbons with 5 H on it so it must be CH₃-CH₂-CH₂. This is confirmed by the presence of a **CH₃ group** Signal 4 (δ=0.95, Integration = 3H, splitting pattern = triplet) thus the CH₂- also has as a neighbour, a carbon with 3 H on it so it must be -CH₂-CH₃.

(Note: 'Atomic Accounting' is now: CH₃ - CH₃ = 0, all atoms accounted for!).



In summary unknown #360-14-B contains an alcohol, and an ethyl group. Therefore, the **Unknown #360-14-B** must be **1-propanol**, which also has the correct/matching molecular weight/formula, and boiling point. Also, the ¹H-NMR Spectrum shows there are 4 different types of H environments present in the molecule (and all 4 have been accounted for).

(Note: The calculation for the 'degrees of unsaturation' in this problem also makes sense. There are 0 degrees of unsaturation present in the unknown.)



Instructor-Led Group ^1H -NMR Analysis Problems

Use the tables below to record your results of the ' ^1H -NMR Spectral Analyses' for the following compounds. Remember to label the signals on the spectrum and assign signal codes to each 'H' in the structure.

Problem 1 $\text{C}_2\text{H}_5\text{I}$ Degrees of Unsaturation = _____


Structure and name	<p style="text-align: center;">ppm</p> <p>HSP-04-351</p>	Structure Assembly Work Area
--------------------	--	------------------------------

^1H -NMR Spectral Data for Instructor-Led Problem 1:
0.5 mL CDCl_3

90 MHz (0.04 mL :


Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
	3.2	2				
	1.8	3				

Problem 2



Degrees of Unsaturation = _____

Structure and name



Structure Assembly Work Area

1H -NMR Spectral Data for Instructor-Led Problem 2:
0.5 mL $CDCl_3$

300 MHz (0.05 g :


Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment

Problem 3



Degrees of Unsaturation = _____

Structure and name




Structure Assembly Work Area

1H -NMR Spectral Data for Instructor-Led Problem 3:

90 MHz (0.02 mL : 0.5 mL $CDCl_3$)


Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment

Problem 4



Degrees of Unsaturation = _____

Structure and name



Structure Assembly Work Area

¹H-NMR Spectral Data for Instructor-Led Problem 4: 90 MHz (0.04 mL : 0.5 mL CDCl₃)

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment

¹H-NMR Analysis Practice Problems

Use the tables below to record your results of the ¹H-NMR Spectral Analyses' of the provided known spectra on this page of the lab manual.

¹H-NMR Spectral Data for:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment

¹H-NMR Spectral Data for:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment

¹H-NMR Spectral Data for:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment

¹H-NMR Spectral Data for:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment

Write-up


Use the following example as a guide for reporting your answer for each of the unknowns.

Experiment 14 Student Name: Phil Terpaper AU ID#: 9876543

Unknown #: 360-14-00 Mol. Wt.: 150.18 Formula: C₉H₁₀O₂ MP or BP: 112-113° C

Number of Degrees of Unsaturation = $nC + 1 + 1/2nN - 1/2 nH - 1/2 nX$
for Compound 360-14-00, C₉H₁₀O₂ = (9) + 1 + 1/2(0) - 1/2 (10) - 1/2(0)
= 9 + 1 - 5 = 5 degrees of unsaturation

Structure: (neatly draw the structure and be sure to indicate and label all the hydrogens)



Name of Compound: *p*-ethylbenzoic acid

Table 1. Infrared Spectral Data:

p-ethylbenzoic acid

<i>p</i> -ethylbenzoic acid	Absorption Band#	Frequency (cm ⁻¹)	Peak Shape (sharp, broad)	Peak Intensity (strong, med. or weak)	Functional Group Indicated
> 3000 cm ⁻¹	A B	2500-3300 3050	v.broad sharp	medium weak	broad carboxylic OH sp ² C-H stretch
Between 3000 and 2000 cm ⁻¹	C	2950	sharp	weak	sp ³ C-H stretch
Between 2000 and 1400 cm ⁻¹	D E	1725 1600,1580,1500	sharp sharp	strong med	C = O carboxylic acid C=C aromatic
< 1400 cm ⁻¹	F	840	sharp	med-strong	<i>para</i> subst. benzene

Functional Group(s) absent: C≡N, C≡C, N-H

Table 2. ¹H-NMR Spectral Data:

Signal #	Shift	Integrat'n	Splitting	Comment	#Neighbour H	Signal Assignment
1	δ 12.1	1H	singlet	(exchanges with D ₂ O)	0	O-H of carboxylic acid
2	δ 8.0	2H	doublet	(Signals 2 and 3 could be referred to together as a multiplet, i.e., 4H)	1	Aromatic H, <i>para</i> substitution pattern
3	δ 7.3	2H	doublet			
4	δ 2.7	2H	quartet		3	Ar-CH ₂ - (shifted downfield)
5	δ 1.3	3H	triplet		2	CH ₃ -C (methyl group)


Discussion and Conclusions:

It was immediately apparent that Unknown # 360-14-00 was a carboxylic acid, -COOH (absorption bands A and D), and an extremely downfield (deshielded) proton (Signal 1 at δ 12.1). The next functional group indicated was an aromatic ring (absorption bands E and B), and downfield protons (Signals 3 and 4 at δ 7.3 and 8.0). Finally, the ¹H-NMR spectrum shows a methyl group (Signal 1). This -CH₃- has as neighbouring carbon with 2 H on it, so it must be CH₃-CH₂. This is confirmed by the presence of a -CH₂ group (Signal 2). The last

piece of the puzzle is that the CH₂ group is shifted downfield ($\delta=2.7$), so it is bonded to an aromatic ring, i.e. an Ar-CH₂-CH₃ group. In summary unknown #360-14-20 contains a carboxylic acid, an aromatic ring and an ethyl group. The substitution pattern appears to be *para*. Therefore, the **Unknown #360-14-00** is ***p*-ethylbenzoic acid**, which also has the correct/matching molecular weight/formula, and boiling point.


Extra Detailed Answer for this Example:

After analyzing and recording the infrared and ¹H-NMR data in Tables 1 and 2 above and consulting the table of 'Chemical Shifts' of the CHEM360 Lab Manual, the following can be deduced. It was immediately apparent that unknown 360-14-00 was a **carboxylic acid**, -COOH, (due to broad infrared absorption bands for C=O and O-H at 1725 and 3300-2500 cm⁻¹ respectively), and an extremely downfield proton (δ 12.1). (Note: 'Atomic Accounting' is now: C₉H₁₀O₂ - CO₂H = C₈H₉).




The next functional group apparent in the unknown was an **aromatic ring** (infrared absorption bands for C=C and sp² C-H at 1600/1500 and 3050 cm⁻¹ respectively), and downfield (deshielded) protons (δ 7.3, 8.0). It is also known that the aromatic ring is di-substituted, in a *para* orientation, hence C₆H₄. This is due to the presence of a single infrared absorption at 840 cm⁻¹ and the pattern of aromatic H δ at 7.3 and 8.0 (displays anisotropy).

(Note: 'Atomic Accounting' is now: C₈H₉ - C₆H₄ = C₂H₅, an 'ethyl group' perhaps?!).




Finally, the spectrum shows a **methyl group** ($\delta=1.3$, Integration = 3H, splitting pattern = triplet). Also this -CH₃- has as neighbouring carbon with 2 H on it so it must be CH₃-CH₂. This is confirmed by the presence of a **CH₂ group** ($\delta=2.7$, Integration = 2H, splitting pattern = quartet) thus the CH₂- also has as a neighbour, a carbon with 3 H on it so it must be -CH₂-CH₃. The last piece of the puzzle is that the CH₂ group is shifted downfield ($\delta=2.7$) and so must be bonded to an aromatic ring, i.e. an **Ar-CH₂-CH₃ group**. (Note: 'Atomic Accounting' is now: C₂H₅ - C₂H₅ = 0, all atoms accounted for!).



In summary unknown #360-14-00 contains a carboxylic acid, an aromatic ring and an ethyl group. The substitution pattern appears to be *para*. Therefore, the **Unknown #360-14-00** must be ***p*-ethylbenzoic acid**, which also has the correct/matching molecular weight/formula, and boiling point. Also, in the ¹H-NMR Spectrum, there are 5 different types of H environments present in the molecule (and all 5 have been accounted for).






Note: The calculation for the 'degrees of unsaturation' in this problem also makes sense (aromatic ring (4) + carbonyl (1) = 5 degrees of unsaturation).

Infrared Spectrum of *p*-ethylbenzoic acid

HIT-NO=6821	SCORE= ()	SDBS-NO=13185	IR-NIDA-25826 : KBR DISC
P-ETHYLBENZOIC ACID			
C ₉ H ₁₀ O ₂			




2977	36	2668	46	1426	21	1064	60	793	41
2956	29	2557	43	1322	18	1020	47	775	58
2940	37	1682	4	1297	10	963	39	747	47
2908	41	1610	21	1236	47	963	38	736	38
2878	39	1572	41	1184	25	858	26	605	50
2836	42	1515	64	1124	36	840	62	547	42
2723	66	1466	66	1061	66	826	60	622	64

(C) (C) (A)

CCc1ccc(cc1)C(=O)O

(B) H (D)

¹H-NMR Spectrum of *p*-ethylbenzoic acid (90 MHz in CDCl₃)



References

1. J.W. Lehman. 1999. Operational Organic Chemistry, 3rd ed., Prentice Hall, Upper Saddle River, NJ.
2. N.A.J. Luff. 1972. DMS Working Atlas of Infrared Spectroscopy, Verlag Chemie, Butterworths & Co. Ltd., London.
3. J. McMurry. 2000. Organic Chemistry, 5th ed., Brooks/Cole Publishing Co., Pacific Grove, California.
4. J. McMurry. 1992. Organic Chemistry, 3rd ed., Brooks/Cole Publishing Co., Pacific Grove, California.
5. Charles J. Pouchert. 1975. The Aldrich Library of IR Data 2nd ed., The Aldrich Chemical Company, Inc., Milwaukee, Wisconsin.
6. C.J. Pouchert and J.R. Campbell. 1974. The Aldrich Library of NMR Spectra, vols 1-11, Aldrich Chemical Company, Inc., Milwaukee, Wisconsin.
7. R.C. Weast et al. 1984. CRC Handbook of Chemistry and Physics 65th ed., CRC Press, Inc., Boca Raton, Florida.

Athabasca University wishes to thank Drs. K. Tanabe and T. Tamura and for all the IR/NMR Spectra used in this manual, obtained from the SDBS web site: <http://www.aist.go.jp/RIODB/SDBS/> (04-Dec-1999).

Answers to Instructor led Group Problems:

Problem 1 iodoethane, 90 MHz (0.04 mL : 0.5 mL CDCl₃) Problem 2 t-butyl alcohol, 300 MHz (0.05 g : 0.5 mL CDCl₃)

Problem 3 propanal, 90 MHz (0.02 mL : 0.5 mL CDCl₃) Problem 4 p-ethoxytoluene, 90 MHz (0.04 mL : 0.5 mL CDCl₃)